

E 980
 Adapter
 for the optimum loading of long hollow instruments.



- For loading long hollow instruments
- Possible to screw in various nozzles or Luer lock adapters

EAN: 4002514312071 / Material number: 07025360 / Old Material Number: 69598001D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, gross height in mm	15
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	80
Net weight in kg	0,04
Gross weight in kg	0,05