


E 705/2

Mobile unit

for the optimum loading of MIS instruments up to 600 mm length.



- Up to 42 hollow instruments can be held
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)
- The vertical clearance for the 3 levels is different

EAN: 4002514313399 / Material number: 07029810 / Old Material Number: 69570503D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7826	•
<b>Product category</b>	
Mobile unit with connection for hot-air drying	•
Mobile unit for MIS instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of MIS instruments	•
<b>Capacity</b>	
Lumened instruments [number]	42
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	530
External dimensions, net width in mm	670
External dimensions, net depth in mm	600
External dimensions, gross height in mm	550
External dimensions, gross width in mm	640
External dimensions, gross depth in mm	650
Net weight in kg	12,50
Gross weight in kg	14,60
<b>Companion pack</b>	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	15
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	5
E 444 drum for fibre optic cables/suction hoses [number]	1
E 445 – 12 covers for E 442 irrigation sleeve, opening with Ø 6 mm [number]	1
E 446 – 12 covers for E 443 irrigation sleeve, opening with Ø 10 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	6
E 449 male adapter, for female Luer Lock [number]	5
E 451 insert 1/6 mesh tray [number]	1
E 452 injector nozzle 2.5 x 60 mm [number]	3
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	8
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	6

E 705/2

Mobile unit

for the optimum loading of MIS instruments up to 600 mm length.



EAN: 4002514313399 / Material number: 07029810 / Old Material Number: 69570503D

E 456 spring stay for MIS instruments [number]	4
E 464 holder for injector nozzle E 454 [number]	3
E 472 spring clip for injector nozzles Ø 4 mm [number]	2
2 m silicone hose, Ø 5 mm	•
Magnetic strip for automatic mobile unit recognition (without magnets) [number]	1