

E 478/1

Holder

for the optimum cleaning of 4 narrow-lumined cannulae (Sautter cannulae).



- Equipped with a stainless steel filter plate (FP-ED)

EAN: 4002514496948 / Material number: 07292920 / Old Material Number: 69747802D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Holder	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Capacity	
Lumened instruments [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0,19
Gross weight in kg	0,23