

ALWD BS 27/28

Scanner

For the simple recording of bar codes for loading and user recognition.



EAN: 4002514753416 / Material number: 07686510 / Old Material Number: 69247012D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Other accessories	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	305
External dimensions, gross width in mm	300
External dimensions, gross depth in mm	200
Net weight in kg	0,90
Gross weight in kg	1,00

- Scanner connection for scanner
- Convenient recording of data for process documentation
- Suitable for the clean, unclean or both sides*

* The combination of scanner and printer is only possible on the clean side.