

E 935/2

Mobile unit

for the optimum cleaning of instruments on 3 levels.



- Loading capacity: 9 DIN mesh trays\*
- Equipped with 2 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets)

\* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002514816432 / Material number: 07765780 / Old Material Number: 69593503D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	524
External dimensions, net width in mm	640
External dimensions, net depth in mm	790
External dimensions, gross height in mm	720
External dimensions, gross width in mm	830
External dimensions, gross depth in mm	780
Net weight in kg	26,50
Gross weight in kg	30,00
<b>Companion pack</b>	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1