

E 417/1

Insert

For the optimum loading of up to 30 ear and nose specula.

EAN: 4002514984605 / Material number: 09033260 / Old Material Number: 69541702D



- suitable for approx. 30 ear and nose specula
- 280 compartments approx. 10 x 10
- The mesh width on the base is 4 mm
- For use in upper and lower basket

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert for ENT instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	280
<b>Area of application</b>	
Reprocessing of ENT instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	80
External dimensions, net width in mm	220
External dimensions, net depth in mm	435
External dimensions, gross height in mm	85
External dimensions, gross width in mm	225
External dimensions, gross depth in mm	465
Net weight in kg	1,30
Gross weight in kg	1,72