

LP 2800

Demineralisation cartridge, empty

For the manual filling of single use resins.



- Max. through flow volume: 800 l/h
- Max water inlet temperature 40 °C – pressure from 0.5-10 bar
- Resin content 19 litres

EAN: 4002515411728 / Material number: 09782590 / Old Material Number: 69431805D

Product classification	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•

LP 2800

Demineralisation cartridge, empty

For the manual filling of single use resins.



EAN: 4002515411728 / Material number: 09782590 / Old Material Number: 69431805D

PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
PG 8581	•
PG 8582	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Water treatment	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	570
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	255
External dimensions, gross width in mm	225
External dimensions, gross depth in mm	600
Diameter in mm	230
Net weight in kg	5,23
Gross weight in kg	7,00