

A 307
Insert
for the optimum loading of 20 theatre shoes up to size 48.



- Suitable for 20 theatre shoes
- Each with 265 mm tall supports
- For use in the lower basket carrier A 151
- Reprocessing of theatre shoes to size 48 possible

EAN: 4002515460856 / Material number: 09862590 / Old Material Number: 69630700D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Capacity	
Theatre shoes [number]	20
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	279
External dimensions, net width in mm	390
External dimensions, net depth in mm	360
External dimensions, gross height in mm	300
External dimensions, gross width in mm	395
External dimensions, gross depth in mm	480
Net weight in kg	1,46
Gross weight in kg	1,99