

A 308

Insert

For the optimum loading of 34 insoles for theatre shoes.



- Suitable for up to 34 insoles
- For use in upper or lower basket carriers

EAN: 4002515460863 / Material number: 09862600 / Old Material Number: 69630800D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0,00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	71
External dimensions, net width in mm	475
External dimensions, net depth in mm	445
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	1,17
Gross weight in kg	2,00