

A 800

Tubular filter

For the safe and economic reprocessing of hollow instruments.



- Filter in ply mesh stainless steel mesh for mobile units
- Mesh size 63 µm
- Reprocessable with cleaning brush A 804

EAN: 4002515538081 / Material number: 09960680 / Old Material Number: 68680000D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Filter	•
<b>Product properties</b>	
Material	Stainless steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	23
External dimensions, net width in mm	23
External dimensions, net depth in mm	405
External dimensions, gross height in mm	45
External dimensions, gross width in mm	220
External dimensions, gross depth in mm	445
Diameter in mm	23
Net weight in kg	0,36
Gross weight in kg	0,88