

PG 8595

Aqua purificator

For two full or partial demineralisation cartridges.



- Suitable for VE/TE P 2800 & VE/TE P 2000
- Suitable for VE/TE P 2800 & VE/TE P 2000
- Possible to hold conductivity meter/ flow meter
- Possible to hold conductivity meter/ flow meter
- Ext. dimen.: H 835 (820 w/o lid), W 300, D 600 mm
- External dimensions: H 835 (820 w/o lid), W 300, D 600 mm

EAN: 4002515558096 / Material number: 10087310 / Old Material Number: 69859501D

Product category	
Fresh-water dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•

PG 8595

Aqua purificator

For two full or partial demineralisation cartridges.



EAN: 4002515558096 / Material number: 10087310 / Old Material Number: 69859501D

PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Water treatment	•
<b>Product properties</b>	
Material	Galvanised plate, Lotus white
Colour	Lotus White
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	835
External dimensions, net width in mm	300
External dimensions, net depth in mm	600
External dimensions, gross height in mm	900
External dimensions, gross width in mm	400
External dimensions, gross depth in mm	740
Net weight in kg	26,50
Gross weight in kg	30,00