


UG 30-90/60-85

Box plinth

to raise the machine incl. side unit, 60 cm deep, 90 cm wide.



- For use in Aqua Purificator PG 8595 and DOS cabinet PG 8596

EAN: 4002515558348 / Material number: 10087360 / Old Material Number: 69309001D

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8504	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Plinth/base	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	300
External dimensions, net width in mm	897
External dimensions, net depth in mm	542
External dimensions, gross height in mm	330
External dimensions, gross width in mm	915
External dimensions, gross depth in mm	563
Net weight in kg	13,40
Gross weight in kg	15,00