

DGG 1/2 - 40 L

Unperforated steam cooking container

For DG 7000 steam ovens and DGM 6000/7000 steam ovens with microwave .

EAN: 4002515603987 / Material number: 10168200 / Old Material Number: 23996120EU1



- For cooking food in sauce, stock or water
- e.g. pasta, rice, soup and casseroles
- Capacity: 2.2 l, usable capacity: 1.6 l
- Exterior dimensions W x D x H: 375 x 197 x 40 mm
- Not suitable for microwave operation!