


UG 30-90/70-85

Box plinth

to raise the machine, 70 cm deep, 90 cm wide.



- Suitable for all machines, width 90, depth 70 cm

EAN: 4002515622506 / Material number: 10210600 / Old Material Number: 69309002D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8582 CD	•
PG 8583 CD	•
<b>Product category</b>	
Plinth/base	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	317
External dimensions, net width in mm	900
External dimensions, net depth in mm	700
External dimensions, gross height in mm	310
External dimensions, gross width in mm	610
External dimensions, gross depth in mm	610
Net weight in kg	15,60
Gross weight in kg	17,00