

ProCare Lab 10 AT - 5 l

Liquid detergent, alkaline, 5 l

For optimum reprocessing of laboratory utensils.



- No phosphorous compounds & oxidants, contains surfactants
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 12.2-12.5

EAN: 4002515660720 / Material number: 10265650 / Old Material Number: 62925101EU1

Product classification	
Washer-disinfectors, laboratory	•
Applications	
Removal of grease	Good
Removal of oil	Good
Removal of resin	Good
Removal of organic dyes	Good
Removal of paraffin	Good
Dispensing recommendation	
Dispensing recommendation per wash cycle	2-4 ml/l (0,2-0,4%)
Intensity of soiling	Strong
pH value of solution	12,2-12,5
Main applications	
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid
Composition	Alkali hydroxide, Polycarboxylate, Tensides
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Not suitable for anodised aluminium	•
Does not contain oxidants	•
Does not contain phosphates	•
Does not contain silicates	•
Not suitable for non-anodised aluminium	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	250
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	150
Gross weight in kg	5
Fill volume in ml	5000
Compliance, certs/symbols	
GHS/CLP compliant	•