

ProCare Lab 10 AT - 10 l

Liquid detergent, alkaline, 10 l

For optimum reprocessing of laboratory utensils.



- No phosphorous compounds & oxidants, contains surfactants
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 12.2-12.5

EAN: 4002515660737 / Material number: 10265660 / Old Material Number: 62925201EU1

<b>Product classification</b>	
Washer-disinfectors, laboratory	•
<b>Applications</b>	
Removal of grease	Good
Removal of oil	Good
Removal of resin	Good
Removal of organic dyes	Good
Removal of paraffin	Good
<b>Dispensing recommendation</b>	
pH value of solution	12,2-12,5
<b>Main applications</b>	
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
<b>Product properties</b>	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Composition	Alkali hydroxide, Polycarboxylate, Tensides
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for laboratory glassware	•
Not suitable for anodised aluminium	•
Does not contain oxidants	•
Does not contain phosphates	•
Does not contain silicates	•
Not suitable for non-anodised aluminium	•
Not suitable for light and non-ferrous metals	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	308
External dimensions, net width in mm	191
External dimensions, net depth in mm	234
External dimensions, gross height in mm	320
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	190
Gross weight in kg	13
Fill volume in ml	10000
<b>Compliance, certs/symbols</b>	
GHS/CLP compliant	•