

ProCare Lab 30 P - 5 l

Neutralising agent, acidic, 5 l

For optimum neutralisation with phosphoric acid based neutralising agent.

EAN: 4002515660454 / Material number: 10266490 / Old Material Number: 62928101EU1



- Free from surfactants and oxidants
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 2.2-1.9

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Applications</b>	
Removal of limescale and other acid-soluble impurities	Good
<b>Dispensing recommendation</b>	
Dispensing recommendation per wash cycle	1-4 ml/l (0,1-0,4%)
pH value of solution	1,9-2,2
<b>Main applications</b>	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
<b>Product properties</b>	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Composition	Phosphate, Phosphoric acid
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Not suitable for anodised aluminium	•
Does not contain surfactants	•
Acidic pre-cleaning	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for non-anodised aluminium	•
Not suitable for light and non-ferrous metals	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	200
External dimensions, gross width in mm	140
External dimensions, gross depth in mm	240
Gross weight in kg	8
Fill volume in ml	5000
<b>Compliance, test certificates, and symbols</b>	
GHS/CLP compliant	•