

ProCare Lab 30 P - 12 kg

Neutralising agent, acidic, 12 kg

For optimum neutralisation with phosphoric acid based neutralising agent.

EAN: 4002515660461 / Material number: 10266500 / Old Material Number: 62928201EU1



- Free from surfactants and oxidants
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 2.2-1.9

Product category	
Washer-disinfectors, laboratory	•
Applications	
Removal of limescale and other acid-soluble impurities	Good
Dispensing recommendation	
Dispensing recommendation per wash cycle	1-4 ml/l (0,1-0,4%)
pH value of solution	1,9-2,2
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Composition	Phosphate, Phosphoric acid
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Not suitable for anodised aluminium	•
Does not contain surfactants	•
Acidic pre-cleaning	•
Not suitable for chrome-plated parts	•
Not suitable for nickel-plated parts	•
Not suitable for non-anodised aluminium	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in mm	308
External dimensions, net width in mm	191
External dimensions, net depth in mm	234
External dimensions, gross height in mm	305
External dimensions, gross width in mm	200
External dimensions, gross depth in mm	245
Weight in g	12000
Gross weight in kg	13
Compliance, test certificates, and symbols	
GHS/CLP compliant	•