

ProCare Lab 11 AP - 10 kg  
 Powder cleaning agent, alkaline, 10 kg  
 For optimum reprocessing of laboratory utensils.

EAN: 4002515650868 / Material number: 10266540 / Old Material Number: 62922101EU1



- Free from phosphorous compounds, oxidising effect
- Specially developed for Miele Professional
- Suitable for all types of materials
- PH value 11.5-11.7

Product classification	
Washer-disinfectors, laboratory	•
Applications	
Removal of inorganic and organic soiling	Good
Dispensing recommendation	
pH value of solution	11,5-11,7
Recommended water range in mmol/l	0-1,1
Main applications	
Water laboratories	•
Serological institutes	•
Phosphate industry	•
Product properties	
Medium type	Alkaline cleaning agent
Aggregate form	Powder
Composition	Oxygen bleach, Tensides
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for porcelain/ceramics	•
Suitable for laboratory glassware	•
Not suitable for anodised aluminium	•
Does not contain phosphates	•
Not suitable for non-anodised aluminium	•
Not suitable for light and non-ferrous metals	•
Dimensions and weight	
External dimensions, net height in mm	268
External dimensions, net width in mm	265
External dimensions, net depth in mm	265
External dimensions, gross height in mm	280
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	270
Weight in g	10000
Gross weight in kg	10
Compliance, certs/symbols	
GHS/CLP compliant	•