

ProCare Lab 11 MA - 10 kg

Powder cleaning agent, mildly alkaline, 10 kg

For optimum reprocessing of laboratory utensils.



- Foam damping effect, free from oxidants
- Suitable for the reprocessing phases in washer-disinfectors
- Suitable for all types of materials
- PH value 11.7-11.9

EAN: 4002515650851 / Material number: 10266720 / Old Material Number: 62921101EU1

Product classification	
Washer-disinfectors, laboratory	•
Applications	
Removal of blood	Good
Removal of protein	Good
Removal of organic residues	Good
Removal of food residues	Good
Dispensing recommendation	
pH value of solution	11,7-11,9
Main applications	
Medical laboratories	•
Biological laboratories	•
Chemical laboratories	•
Product properties	
Medium type	Mildly alkaline Cleaning agent
Aggregate form	Powder
Composition	
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Limited suitability for non-ferrous metals	•
Not suitable for anodised aluminium	•
Reduces foam	•
Does not contain oxidants	•
Limited suitability for non-anodised aluminium	•
Dimensions and weight	
External dimensions, net height in mm	268
External dimensions, net width in mm	265
External dimensions, net depth in mm	265
External dimensions, gross height in mm	280
External dimensions, gross width in mm	270
External dimensions, gross depth in mm	270
Weight in g	10000
Gross weight in kg	10
Compliance, certs/symbols	
GHS/CLP compliant	•