

A 810

Cover net

For optimum securing of delicate items.



- Stainl. steel frames w/ plastic netting, dimensions 445 x 215 mm
- Particularly resistant and durable

EAN: 4002515676783 / Material number: 10283650 / Old Material Number: 69501700D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Cover	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, White
Area of application	
Reprocessing of surgical instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	20
External dimensions, gross width in mm	235
External dimensions, gross depth in mm	470
Net weight in kg	0,44
Gross weight in kg	0,57