


A 811/1  
Replacement net  
For replacing the net on an A 810/1



- Easy replacement of the netting
- Particularly resistant and durable

EAN: 4002515676493 / Material number: 10283690 / Old Material Number: 69502000D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Cover	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, White
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	20
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	230
Net weight in kg	0,02
Gross weight in kg	0,07