

E 454/1

Injector nozzle

for the optimum cleaning of diverse loads, particularly for proctoscopes.



- Suitable for proctoscopes. Trocar sleeves and other items

EAN: 4002515722954 / Material number: 10342490 / Old Material Number: 69745402D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7824	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Injector nozzle	•
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of lumened instruments	•
Capacity	
Lumened instruments [number]	1
Dimensions and weight	
External dimensions, net height in mm	8
External dimensions, net width in mm	8
External dimensions, net depth in mm	300
External dimensions, gross height in mm	60
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	350
Diameter in mm	8
Net weight in kg	0,06
Gross weight in kg	0,07