


A 862
Holder rack
for injector modules A 302/2 and A 302/3



- Arrangement: 4 x 8 fields
- For use on A 302/2 and A 302/3 injector modules
- For the optimum positioning of laboratory glassware
- Plastic-coated

EAN: 4002516087687 / Material number: 11056820 / Old Material Number: 69686200D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black
Number of compartments	32
Compartment width in mm	56
Compartment length in mm	52
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0,00-0,00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	30
External dimensions, net width in mm	250
External dimensions, net depth in mm	350
External dimensions, gross height in mm	35
External dimensions, gross width in mm	300
External dimensions, gross depth in mm	400
Net weight in kg	0,10
Gross weight in kg	0,18