

A 840/1

Injector nozzle

with plastic supports, nozzle diameter 6 mm, length 130 mm, 5 items.



- Suitable for mobile injector units & baskets screw connection
- Part of the EasyLoad system
- Packaging unit with 5 nozzles
- Ideal for load items with neck diameters of between 12 mm and 85 mm
- Diameter of cap at top: 11 mm

EAN: 4002516188803 / Material number: 11112600 / Old Material Number: 69684001D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	130
External dimensions, net width in mm	92
External dimensions, net depth in mm	86
External dimensions, gross height in mm	107
External dimensions, gross width in mm	130
External dimensions, gross depth in mm	195
Net weight in kg	0,13
Gross weight in kg	0,24