

A 841/1

Injector nozzle

with plastic supports, nozzle diameter 6 mm, length 210 mm, 5 items



- Suitable for mobile injector units & baskets screw connection
- Part of the EasyLoad system
- Packaging unit with 5 nozzles

EAN: 4002516188810 / Material number: 11113130 / Old Material Number: 69684101D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	210
External dimensions, net width in mm	92
External dimensions, net depth in mm	86
External dimensions, gross height in mm	108
External dimensions, gross width in mm	117
External dimensions, gross depth in mm	305
Net weight in kg	0,17
Gross weight in kg	0,27