

RX2/RX3-YCR-MT12  
RX2/RX3 120 lithium-ion battery  
for non-stop 120 min. running time of the Scout RX2/RX3.

EAN: 4002516435044 / Material number: 11779170 / Old Material Number: 4199701EU1



- High-quality, high-capacity battery:  
75.5 Wh
- Long service life due to Li  
technology without memory effect
- Protection against short-circuits,  
overheating & overvoltage
- Replacing the battery is quick and  
easy
- Every high-performance cell is  
tested separately