

WA SO 1503 L
UltraSoft fabric conditioner 1.5 l
for soft, bouncy laundry and a fresh Aqua fragrance.

EAN: 4002516531159 / Material number: 11986980 / Old Material Number: 11997233EU2



- Natural fragrance of pure cleanliness and freshness
- Minimises creasing and makes ironing laundry easier
- Bottles made from 100 % recycled plastic – eco-friendly
- Free from microplastics and biodegradable
- Prevents electrostatic charges in washer dryers
- Economical – for 50 laundry loads