

WA UP2 1403 L

Miele UltraPhase 2

universal UltraPhase 2 suitable for every UltraPhase 1.

EAN: 4002516779681 / Material number: 12517030 /
Old Material Number: 11997298EU2



- Active oxygen for gentle and effective stain removal
- OxiCare for thorough care and cleaning of fabrics and appliances
- WhiteBooster for radiant white laundry
- Up to 30 % less consumption thanks to automatic dispensing
- Very economical – 1.4 litres for 50 wash loads
- For use in all Miele TwinDos washing machines and washer-dryers