

## Guard L1

### Cylinder vacuum cleaner

powerful | premium design | convenient cleaning with 2-in-1 nozzle



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years\*\* use
- Ideal power for every cleaning task: 4 selectable modes
- Individually selectable height and ergonomic Comfort handle
- Premium design concept with stylish applications
- Range extra, fewer socket changes: 12 m operating radius
- 2-in-1 upholstery and crevice nozzle for easy cleaning without switching over
- AirClean filter: effective filtration for clean room air

\* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516862727 / Material number: 12734280 /  
Old Material Number: 41UDF020ZA

<b>Construction type and design</b>	
Appliance colour	Nordic blue
<b>Construction type</b>	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard L1
Motor type	PowerLine
Product type	SUDFO
<b>User convenience</b>	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Suction tube	Suction tube X(HES)
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
Comfort handle 2.0	•
<b>Mobility</b>	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	7,5
<b>Filter system</b>	
Vacuum cleaner bag type	HyClean Pure TU
Dustbag volume in l	3,5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
<b>Safety</b>	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
<b>Technical data</b>	
Power	890
Volume in dB	74,0
Width in mm	500
Height in mm	239
Weight in kg	6,6
<b>Accessories included</b>	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead SBD 365-3	•
Vacuum cleaner bag voucher	•