

Guard L1 Cat & Dog
Cylinder vacuum cleaner
powerful | attachment for removing pet hair | ideal for pet owners



- Active AirClean filter neutralises unpleasant odours
- Durable: tested for the equivalent of 20 years** use
- ComfortFit: developed for the highest hygiene standards
- Individually selectable height and ergonomic Comfort handle
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal power for every cleaning task: 4 selectable modes
- Premium design concept with stylish applications
- Turbobrush: removes stubborn dirt and animal hair

* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516862734 / Material number: 12734290 /
Old Material Number: 41UZFO34ZA

Construction type and design	
Appliance colour	Obsidian black
Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard L1
Motor type	PowerLine
Product type	SUZFO
User convenience	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Suction tube	Suction tube X(HES)
Park system for pauses during vacuuming	•
Comfort cable rewind	•
Comfort handle 2.0	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	7,5
Filter system	
Vacuum cleaner bag type	HyClean Pure TU
Dustbag volume in l	3,5
Dustbag change indicator	•
Exhaust filter	Active AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
Technical data	
Power	890
Volume in dB	74,0
Width in mm	500
Height in mm	239
Weight in kg	6,6
Accessories included	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead SBD 365-3	•
Turbobrush STB 305-3 (DOUBLE; LOCKED)	•
Vacuum cleaner bag voucher	•