


## KWT 6722 iGS-1

Ugradbeni hladnjak za temperiranje vina  
s FlexiFrame, Sommelier setom i Push2open za najzahtjevnije ljubitelje vina.



EAN: 4002516775799 / Materijal: 96008140 / Stari broj materijala: AT96008140

### Upravljanje


Nezavisna regulacija temperature zone za hlađenje i zone za zamrzavanje	Ne
---	----

- Fleksibilno pohranjivanje većih boca za vino uz FlexiFrame
- Savršen dizajn za kuhinje bez ručki – Push2open
- Paralelno pohranjivanje vina odvojenim Temperaturnim zonama
- Nema mirisa u hladnjaku zahvaljujući Active AirClean filtrima
- Vaš hladnjak za vino savršeno se uklapa – s LED-osvjetljenjem


## KWT 6722 iGS-1

Ugradbeni hladnjak za temperiranje vina


s FlexiFrame, Sommelier setom i Push2open za najzahtjevnije ljubitelje vina.



KWT6722iGS, KWT6722iGS-1, installation draw



KWT6722iGS, KWT6722iGS-1, side view (installation drawing)



KWT6722iGS, KWT6722iGS-1, Interior dimensions (Installation drawings)