

KWT 2612 Vi
Acondicionador de vino MasterCool
para diseño y tecnología superior en formato grande.


Código EAN: 4002516322955 / Número de material: 11503610 /
Número de material antiguo: 36261205USA




- Apertura de la puerta cómoda gracias a Push2open
- Varios vinos perfectamente atemperados: tres zonas de temperatura
- Almacenamiento más flexible de botellas grandes con FlexiFrame
- Interior luminoso gracias a los elementos LED: BrilliantLight
- Comunicación sencilla con el aparato a través de WiFiConn@ct

KWT 2612 Vi


Acondicionador de vino MasterCool para diseño y tecnología superior en formato grande.




Opening angle 90°/115°, right side stop, MasterCool (Sketch) with footnote




Side-by-side, MasterCool, without partition wall, with heating mat (Sketch) (with footnote)




KWT2611Vi, KWT2671ViS, KWT2601Vi, KWT2661ViS, installation, MasterCool, installation drawings




Sideview, MasterCool, Sketch




Side-by-Side, MasterCool, footnote, installation drawings



KFP2455ed/cs, KFP2454ed/cs, front panel, front dimensions, MasterCool, installation drawings



Installation, MasterCool, *footnote, installation drawings



KWT2602Vi/SF, KWT2612Vi/SF, KWT2601Vi/SF, KWT2611Vi/SF, inner space, MasterCool, installation drawings